

MedPark



**NEW BONE  
NEW LIFE**

**DENTAL**  
Clinical case



# CONTENTS



<b>CLINICAL CASE (+BIOPSY)</b>	<b>  1</b>	<b>Alveolar Ridge Preservation without Membrane</b>
<b>CLINICAL CASE (+BIOPSY)</b>	<b>  2</b>	<b>Alveolar Ridge Preservation</b>
<b>CLINICAL CASE (+BIOPSY)</b>	<b>  3</b>	<b>Alveolar Ridge Preservation without Membrane</b>
<b>CLINICAL CASE (+BIOPSY)</b>	<b>  4</b>	<b>Alveolar Ridge Preservation</b>
<b>CLINICAL CASE (+BIOPSY)</b>	<b>  5</b>	<b>Sinus Floor Elevation (Lateral Approach)</b>
<b>CLINICAL CASE (+BIOPSY)</b>	<b>  6</b>	<b>Vertical &amp; Horizontal Bone Augmentation</b>
<b>CLINICAL CASE (+BIOPSY)</b>	<b>  7</b>	<b>Regeneration of the Odontogenic Cyst Area</b>

CLINICAL CASE	8	Horizontal Alveolar bone Augmentation
CLINICAL CASE	9	Minor Bone Augmentation without Membrane
CLINICAL CASE	10	Moldable Augmentaion in Anterior Area
CLINICAL CASE	11	Alveolar Ridge Preservation
CLINICAL CASE	12	Periodontal Defect Management
CLINICAL CASE	13	Socket Management in Posterior Area
CLINICAL CASE	14	Dehiscence Defect in Anterior Area
CLINICAL CASE	15	Open Socket Management in Posterior Area
CLINICAL CASE	16	Narrow Rigde Augmentation in Anterior Area
CLINICAL CASE	17	Easy Augmentation of Dehiscence Defect
CLINICAL CASE	18	Simplified Vertical Augmentation for Advanced Bone Resorption
CLINICAL CASE	19	Dough Staged Sticky Bone Graft
CLINICAL CASE	20	Immediate Implant Placement
CLINICAL CASE	21	Easy 3D Ridge Augmentation
CLINICAL CASE	22	Labial Fenestration Defect
CLINICAL CASE	23	Minimally Invasive Bone Grafting in Anterior Area
CLINICAL CASE	24	Horizontal bone Augmentation
CLINICAL CASE	25	Labial Dehiscence in Upper Anterior Area
CLINICAL CASE	26	Labial Fenestration in Lower Anterior Area
CLINICAL CASE	27	Peri-implantitis Treatment with laser therapy
CLINICAL CASE	28	Extraction Socket Management
CLINICAL CASE	29	Moldable Augmentaion in Posterior Area
CLINICAL CASE	30	Immediate implant placement
CLINICAL CASE	31	Immediate implant placement
CLINICAL CASE	32	Sinus augmentation
CLINICAL CASE	33	Socket Management in Posterior Area & Immediate implant placement
CLINICAL CASE	34	Minor bone augmentation (Maxillary sinus lift revision surgery)
CLINICAL CASE	35	Horizontal Alveolar bone Augmentation
CLINICAL CASE	36	Horizontal ridge augmentation
CLINICAL CASE	37	Horizontal guided bone regeneration
CLINICAL CASE	38	Horizontal Alveolar bone Augmentation
CLINICAL CASE	39	Narrow Rigde Augmentation in Posterior Area

# Alveolar Ridge Preservation without Membrane


**Initial assesment** 78 years old female patient has suffered from gingival swelling & pain due to Vertical root fracture of Lt. 2nd premolar.  
**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› Extensive bone destruction even apical & buccal aspect was found in CT scan. Thus, Alveolar ridge preservation procedure with S1 was required prior to place implant

## Conclusions

› Despite of single use of S1 without membrane covering for ARP, Extra-fine soft tissue healing and full coverage was found in a month. And ridge contour was preserved for 5 month of healing time. histologic evaluation revealed extra-fine new bone formation in quantity and quality of regenerated bone



Preoperative x-ray



Tooth extraction




Application of S1 bone graft material and suture (Open socket)




Post OP (2 weeks)



Post OP (2 months)




Post OP (5 months)



Post OP (5 months), Implant placement



Post OP (5 months), suture




Post OP (5 months), Implant placement



Post OP (7 months), Prothesis



Post OP (7 months), Prothesis




Alveolar bone changes on CBCT

**Biopsy time** 5 month after S1 bone graft material use

**Biopsy method** Collect at 6mm depth by using Trephine bur from #13(universal system) implant placement area

**Findings** Despite of single use of S1 without membrane covering for ARP, histologic evaluation revealed extra-fine new bone formation in quantity and quality of regenerated bone

› Hematoxylin & Eosin staining / Masson's Trichrome staining



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	0.65	8.96
<b>New bone area</b>	<b>3.43</b>	<b>47.24</b>
Others	43.80	43.80
Total bone area	4.08	56.20



# Alveolar Ridge Preservation

**Initial assesment** Advanced Peri-implantitis was found on Implant Br. Of Lt. maxilla with buccal bone defect, 41 years old male patient who suspected heavy smoker.

**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

› After Fixture removal, Severe buccal dehiscence defect & poor soft tissue condition enabled only Alveolar Ridge Preservation on Lt. Premolar area. S1 bone was applied on dehiscence defect and collagen membrane covering was done


## Conclusions

› Not only volumatic change was found in CT scan, but also clinical situation for implant placement was acceptable. Additional bone graft was done simultaneously with Implant placement 5 month after ARP. Fine bone regeneration was found on histologic evaluation


Before



After




Before



After




Before



After



# Treatment Steps




Preoperative x-ray



Peri-implantitis condition




Defect after old fixture remove



Huge Defect after fixture remove



Application of S1 bone graft material



Application of S1 bone graft material



Bone covering with granulation tissue with collagen membrane



Post-OP (5 months)



Post-OP (5 months)




Post OP (5 months), Implant placement



Post OP (5 months), Implant placement



Post OP (9 months), Prothesis



Alveolar bone changes on CBCT




**Biopsy time** 5 month after S1 bone graft material use

**Biopsy method** Collect at 6mm depth by using Trephine bur from #12(universal system) implant placement area

**Findings** Although it was a wide defect, Fine bone regeneration was found on histologic evaluation

#### › Hematoxylin & Eosin staining



#### Total Region of Interest Area

	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	1.47	26.77
<b>New bone area</b>	<b>1.46</b>	<b>26.64</b>
Others	2.55	46.59
Total bone area	2.93	53.41

# Alveolar Ridge Preservation without Membrane

**Initial assesment** 47 years old male patient's got a floating tooth on Rt. Mn. Canine with advanced apical periodontitis. Alveolar Ridge Preservation was required for 2 wall Extensive bony defect and poor soft tissue condition. And more than 6 month of healing period for implant placement was expected



**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› After Tooth extraction, thick Inflamed soft tissue was peeled off carefully from underlying bone. Then, S1 bone application was done on huge defect, And covered with granulation tissue that play a role of protective membrane for graft material


## Conclusions

› For 5 months, Soft tissue and bony contour was well preserved in CT scan. But regenerated bone shows soft bone quality despite of histologic result was acceptable bone regeneration. More than 6 month of healing period would be required for implant placement in case of 1-wall or 2-wall defect







# Treatment Steps




Preoperative x-ray



Huge defect after Tooth extraction



Application of S1 bone graft material



Bone covering with granulation tissue without membrane



Post-OP (2 weeks)



Post-OP (6 weeks)



Post-OP (5 months)




Post-OP (5 months)




Post OP (5 months), Implant placement



Post OP (5 months), Implant placement



Post OP(5months)



Alveolar bone changes on CBCT



Post OP (8 months), prosthesis



Post OP (8 months), prosthesis




Post OP (9 months)

**Biopsy time** 5 month after S1 bone graft material use

**Biopsy method** Collect at 6mm depth by using Trephine bur from #27(universal system) implant placement area

**Findings** In case of 1-wall or 2-wall defect More than 6 month of healing period would be required for implant placement, regenerated bone shows soft bone quality and histologic result was acceptable for implantation

#### › Hematoxylin & Eosin staining



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	1.21	16.00
<b>New bone area</b>	<b>1.52</b>	<b>20.16</b>
Others	4.835	63.83
Total bone area	2.73	36.17



# Alveolar Ridge Preservation

**Initial assesment** Visited for peri-implantitis of the mandibular left first and second molars Removal of existing implant and inflammation curettage through guided bone regeneration surgery. After observing the results, perform implant placement


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Where there is insufficient support space for implant placement due to a creation of a wide socket range after removal of a peripheral inflammatory tissue due to peri-implantitis
- › On the day of GBR surgery, after soft tissue elevation, traces of the bone graft material used for previous implant placement were observed.

## Conclusions

- › At 6 months after S1 bone graft material use, the boundary line between the S1 bone graft material and the natural alveolar bone is naturally connected to generate new bone
- › Biopsy results of the new bone formed in the implant placement site confirmed that osteocytes were generated around the graft material and connected to form a mature bone



1 Preoperative X-ray (#36, #37)



2 Implant removal and inflammation curettage



3 Application of S1 bone graft material



4 Application of S1 bone graft material



5 Intraoral photo at 6 months after application of S1 bone graft material



6 Soft tissue incision for implant placement



7 Implant placement in newly formed bone



8 Suture after connecting healing abutment




9 X-ray after implant placement (6 months after S1 bone graft material use)

**Biopsy time** 6 months after S1 bone graft material use  
**Biopsy method** Collect at 6mm depth by using Trephine bur from #36 implant placement area


**Findings**

- › After 6 months, the bone formation rate was confirmed to be 60.29%, and most forms a network with mature bone.
- › It is confirmed to be safe as a graft material, as there are no signs of inflammation or immune rejection
- › Osteogenesis was formed in a form surrounded by bones around the graft material, and the uniform distribution of osteocyte was observed, confirming the mature stage of ossification.

› Hematoxylin & Eosin staining



› Masson's Trichrome staining



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	1.44	20.99
<b>New bone area</b>	<b>4.13</b>	<b>60.29</b>
Others	1.28	18.72
Total bone area	5.57	81.28

※ Biopsy analysis : OBen tissue analysis



# Sinus Floor Elevation (Lateral Approach)



**Initial assesment** A patient who had been using dentures for a long time. Visited for implant placement. Performed vertical and horizontal bone augmentation at the base of the maxillary sinus by using a lateral approach as there is 1 mm of residual bone in the maxillary left posterior region near the maxillary sinus.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Insufficient condition with 1mm of vertical residual bone in the upper left molar area
- › Approach maxillary sinus by forming a window through a lateral approach
- › S1 bone graft material graft after maxillary sinus elevation

## Conclusions

- › A case that shows the formativeness of a space by lifting the sinus floor membrane using S1 bone graft material
- › Hard bone is formed in which resistance can be felt on Trephine burs when collecting for biopsy of the S1 bone graft material area
- › As a result of the biopsy, a new bone was well formed around the graft material and mature bone tissue could be observed.



Preoperative X-ray



Intraoral photo before left maxilla graft




Application of S1 bone graft material from the buccal side of the remaining bone in left maxilla




X-ray immediately after using S1 bone graft material




Implant placement in #26 area at 4 months after using S1



Biopsy request for posterior region of #26 at 11 months after using S1



Final prosthetics



X-ray after 26 months




After 33 months of using S1 bone graft material

**Biopsy time** 11 months of S1 bone graft material use  
**Biopsy method** Collect at 6 mm depth by using Trepine bur from the posterior side of the implant placement area

- Findings**
- › It is considered safe as a graft material, as there are no signs of inflammation or immune rejection
  - › High-quality bone is formed uniformly from the host bone to the maxillary sinus mucosa
  - › Most of the new bone is connected to the graft material, and the newly formed bone forms a very dense network

› Hematoxylin & Eosin staining



› Comparison of S1 and Product ‘B’ for bone formation rate

	S1	Product ‘B’
Graft area	Sinus	Sinus
Observation period	11 months	12 months
Bone formation rate	<b>28.22 %</b>	26.60 %

※ Biopsy analysis : Knotus, Tissue analysis by Wonkwang University

**Reference**

(1) Son WK, Shin SY, Yang SM, Kye SB. Maxillary sinus floor augmentation with anorganic bovine bone: Histologic evaluation in humans. J Korean Acad Periodontol. 2009;39(1):95-102.

(2) Lee YM, Shin SY, Kim JY et al. Bone reaction to bovine hydroxyapatite for maxillary sinus floor augmentation: Histologic results in humans. Int J Periodontics Restorative Dent. 2006;26:471-481.



# Vertical & Horizontal Bone Augmentation



**Initial assesment** Situation requiring extraction of the entire maxillary tooth and removal of the inflamed tissue due to periodontitis in the oral cavity Because extensive horizontal and vertical bone reconstruction of the maxillary alveolar bone was required, implant placement was planned at 4 months after bone augmentation

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (Non-absorbable) ☐ No

## Objectives

- › The loss of alveolar bone due to periodontitis progressed over a wide range, and large bone defects were observed in the maxillary right canine and the first premolar area.
- › Horizontal and vertical guided bone augmentation using S1 and a non-absorbable membrane.

## Conclusions


- › Using a guided bone regeneration technique that has been widely performed, the biopsy shows good ossification and good bone quality where the boundary between the bone graft and the new bone cannot be observed.
- › When soft tissue was elevated After 4 months of using the bone graft material, there was sufficient bone quality and bone width for implant placement.



Preoperative X-ray



Postoperative X-ray after using S1 bone graft material



CBCT after using S1 bone graft material



Intraoral photo after 3 months (lateral right maxilla)




Intraoral photo after 3 months



Soft tissue incision after 4 months for implant placement



Abutment connection for Final prosthetics



X-ray after Ffinal prosthetics




Final prosthetics

**Biopsy time** 4 months of S1 bone graft material use  
**Biopsy method** Collect at 6mm depth using Trepine bur from #36 implant placement area

- Findings**
- › Ossification was good enough that the boundary between the bone graft material and the new bone was not visible, and bone quality was good
  - › It has been confirmed as a mature ossification, stage in which osteoblasts are observed in the bone graft material, further than the initial stage of bone formation in which osteoblasts are active
  - › It is considered safe as a graft material, as there are no signs of inflammation or immune rejection
  - › The bone formation is very good with densely attached connection is densely attached as the bone graft and the new bone have good affinity

› Hematoxylin & Eosin staining



› Comparison of S1 and Product ‘B’ for bone formation rate

	S1	Product ‘B’
Graft area	Alveolar ridge	Alveolar ridge
Observation period	4 months	6 months
Bone formation rate	<b>25.23 %</b>	18.3 %

※ Biopsy analysis : OBen tissue analysis

**Reference**

- (1) Lee YM, Shin SY, Kim JY et al. Bone reaction to bovine hydroxyapatite for maxillary sinus floor augmentation: Histologic results in humans. Int J Periodontics Restorative Dent. 2006;26:471-481.
- (2) Nicola U. Zitzmann, Dr Med Dent\* Peter Schärer, Prof Dr Med Dent, MS\*\* Carlo P. Marinello, Prof Dr Med Dent, MS\*\*\* Peter Schüpbach, Dr Sc Nat, PhD\*\*\*\* Tord Berglundh, DMD, PhD\*\*\*\*\*
- (3) Zitzmann NU, Schärer P, Marinello CP, Schüpbach P, Berglundh T. Alveolar ridge augmentation with Bio-Oss: a histologic study in humans. Int J Periodontics Restorative Dent. 2001 Jun;21(3):288-95. PMID: 11490406.



# Regeneration of the Odontogenic Cyst Area

**Initial assesment** For reconstruction of a bone defect caused by an inflammatory cyst at the root of the mandibular right first molar tooth, it is planned that guided bone reconstruction surgery by using S1 bone graft material will be performed first, and then implant placement will be performed after confirming bone tissue formation




**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives


- › Implant placement plan for the right mandibular first and second molars.
- › Observation of a large cyst and inflamed tissue below the 1st molar apex.
- › Planned implant placement after extraction and guided bone reconstruction in the affected area.

## Conclusions


- › After 3 months of using S1 bone graft material, a dense density of new bone was observed on radiographs.
- › Biopsy results of the new bone formed in the implant placement site confirmed the new bone induction process
- › Healthy new bone including blood vessels is formed even in the areas with large defects



# Treatment Steps



Preoperative X-ray (before extraction of #46)



X-ray after extraction of mandibular left first molar



Intraoral photo after extraction of mandibular left first molar



Residual bone condition after soft tissue flap elevation




Application of S1 bone graft material (#46)




Application of S1 bone graft material (#47)



X-ray at 3 months after application of S1 bone graft materia



Postoperative CBCT cross-section



Implant placement at 4 months after S1 bone graft material use



Postoperative Intraoral photo after final prosthetics



Final prosthetics (#46, #47)




X-ray after final prosthetics

**Biopsy time** 4 months after S1 bone graft material use


**Biopsy method** Collect at 6mm depth by using Trephine bur from #36 implant placement area

- Findings**
- › 22.3% of the new bone area in the preservation of the socket with large defect, which is a meaningful value
  - › A pattern in which a lot of new bone is formed between the safe graft material and the bone graft material without any special inflammatory reaction
  - › Osseointegration expected to be excellent
  - › A pattern in which very healthy new bone including blood vessels is formed

› **Hematoxylin & Eosin staining**



› **Masson's Trichrome staining**



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	2.10	11.46
<b>New bone area</b>	<b>4.09</b>	<b>22.34</b>
Others	12.13	66.18
Total bone area	6.19	33.81

※ Biopsy analysis : OBen tissue analysis



# Horizontal Alveolar bone Augmentation

**Initial assesment** A 63-year-old woman visited the dental clinic because of uncomfortable dentures and wanted implants.



**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

- › Extraction of residual canine teeth and implantation of implants in the anterior and posterior areas
- › Bone graft due to resorption of labial side of anterior teeth

## Conclusions

- › Graft S1 bone into the thin labial of the anterior teeth, and visually observe the augmented alveolar bone and healed soft tissue after 3 months through CBCT





# Treatment Steps




Preoperative Panorama




Preoperative IntraOral Photo




Preoperative CBCT



Extraction of the remaining upper right canine tooth and partial incision of the anterior tooth



Transplantation of S1 bone graft material



Use of COLLA Membrane



Suture



Preoperative Panorama




Intraoral photograph of soft tissue healing 3 months after surgery



Final Panorama



Final Prosthesis



Postoperative CBCT

# Minor Bone Augmentation without Membrane

**Initial assesment** Edentulous Rt. Mx. premolar area of 67 year old female patient was not sufficient horizontal volume for implant placement. Only bone graft was required but patient did not want to get a Extensive GBR for her DM condition


**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› After implant placement, buccal fenestration defect was found on 1st premolar implant. S1 bone was applied on buccal defect without membrane or additional fixation for space maintenance. Minor flap advancement and primary closure was done with only 5-0 nylon suture

## Conclusions

› Abundant bone regeneration was found after 3 Months, even regenerated bone covered over the cover screws. Regenerated buccal bone volume was stable over 8 months



Preoperative x-ray



Implant placement



Fenestration defect



Suture



Post OP (3 weeks)



Post OP (3 months), 2nd OP



Post OP (3 months), 2nd OP




Post OP (3 months), 2nd OP




Post OP (4 months)



Post OP( 4 months), Final prosthesis



Post OP (4 months), Final prosthesis



Alveolar bone changes on CBCT



# Moldable Augmentation in Anterior Area

**Initial assesment** A 66-year-old male patient visited with complaints of a strong odor around the maxillary anterior prosthesis (4-unit bridge) and wished to have implant treatment. Severe alveolar bone contraction progressed due to root inflammation, and the residual bone buccal width of 3 mm was observed on CBCT.



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › In need for extensive alveolar bone augmentation as buccal width of alveolar bone is 3mm for extraction of residual maxillary teeth (#11, #21, #22) and removal of inflamed granulation tissue.
- › When observed with naked eye, there is a defect in the form of fenestration and dehiscence of the alveolar bone after soft tissue elevation

## Conclusions

- › Observation of good healing of the soft tissue around bone graft material although the surgical area was wide.
- › Upon CBCT, the new bone around the implant is naturally connected to the existing alveolar bone.
- › Although no membrane was used, the position of the bone graft material was maintained stably, and bone augmentation surgery was successfully performed relatively easily.






# Treatment Steps



Preoperative panorama



Preoperative CBCT



Preoperative Intraoral photo



Implant placement at #11, #21, #22



Dehiscence defect and broad bone loss



Application of S1 bone graft material



Application of S1 bone graft material




Suture



6 weeks after S1 graft



X-ray after final prosthetics



CBCT after final prosthetics



Final prosthetics

# Alveolar Ridge Preservation

**Initial assesment** Visited for implant placement in the tooth loss area, and planned vertical bone augmentation by using bone graft material at the same time as implant placement in the mandibular left first major molar defect

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No


## Objectives

- › Defect area on #36 shows the form of vertical and horizontal alveolar resorption.
- › Plan for vertical bone augmentation at the same time as implant placement in order to reduce the difference between the adjacent teeth and excessive bone level.


## Conclusions

- › Vertically augmented alveolar bone around the #36 implant can be confirmed on comparing the preoperative X-ray and postoperative X-ray of the final prosthesis placed 3 months after surgery.
- › Because of the augmented alveolar bone and soft tissue, it was possible to fabricate a final prosthesis in harmony with the adjacent teeth, and it is expected to be advantageous in terms of future periodontal management.

Before



After



Before



After



Before




After





# Treatment Steps



Preoperative X-ray



Preoperative Intraoral photo



Implant placement after soft tissue flap elevation and inflammation curettage



Application of S1 bone graft material on implant bone loss area



Vertical and horizontal bone augmentation by using S1



Postoperative X-ray



10 weeks after the surgery



Observation of new bone formation upon second surgery at 10 weeks after 10 weeks of surgery




Abutment connection for prosthetics



Final Prosthetics delivery at 3 months after S1 use



Final Prosthetics delivery at 3 months after S1 use



Before surgery Right after surgery 4months

Alveolar bone changes on CBCT



# Periodontal Defect Management

**Initial assesment** Advanced periodontal disease was found on Lt. mandible molar with buccal bone defect, 57 years old

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› After removal of the inflammatory granulation tissue, destruction of the alveolar bone connected to the buccal side of the two molars was observed, and there was a large puddle-shaped defect around the mesial root of the second molar. It is planned to use only the S1 bone graft material in the area where it is difficult to use a membrane in the area of the defect.

## Conclusions

› Although the membrane was not used, the volume of the bone graft material was maintained well and the new bone was well formed on the CT scan.  
› The process of new bone formation was observed through radiopacity on CBCT, and after 2 months, it merged with the surrounding bone tissue.




Preoperative x-ray



Flap operation



Inflammation curettage



Defect after removal of inflamed tissue



Ready to use S1 bone graft material



Application of S1 bone graft material



Preoperative x-ray



Suture



Post OP



Post OP (3 months)



Post OP (3 months)



Alveolar bone changes on CBCT

# Socket Management in Posterior Area

**Initial assesment** A patient with severe periodontitis in the maxillary left first molar with little or no bone remaining at the site for implant placement. Performed bone augmentation of the tooth extraction area


**Products** S1 bone graft material (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (PRF) ☐ No

## Objectives

- › Alveolar bone condition after maxillary first molar
- › Confirmation of extraction and defect area after soft tissue elevation
- › S1 was hydrated by collecting PRF, and PRF membrane was applied

## Conclusions


- › Even though the defect area was large, the shape was well maintained without a membrane, and new bone was formed well.
- › Good results were obtained only by using residual bone and S1 graft material without maxillary sinus elevation
- › Bone formation was good in which the height of the marginal bone around the implant harmonizes with the adjacent teeth



Preoperative X-ray




Preoperative Intraoral photo



Hydration of S1 using PRF



Application of S1 bone graft material after implant placement




X-ray at 1 months after S1 bone graft material use




Bone formation upon second surgery at 7 months after placement



Abutment placement before prosthetics



X-ray after final prosthetics



X-ray at 10 months after S1 bone graft material use



After months of using S1 bone graft material



# Dehiscence Defect in Anterior Area

**Initial assesment** Planning to have horizontal bone augmentation surgery with the remaining net lingual width of 2.8mm and implant placement plan, after extraction of the maxillary left and right lateral incisors


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extraction of maxillary left and right lateral incisors due to apical inflammation
- › Horizontal bone augmentation was performed with a lingual width of 2.8 mm in the remaining bone.
- › Application of S1 graft material after implant placement with a diameter of 3.5mm

## Conclusions


- › Though the bone width of the defect part was insufficient and the range was wide, the desired shape was well made after using the S1 bone graft material, and the position of the graft material was also stably maintained.
- › The position of S1 bone graft material was maintained without use of a membrane, triggering new bone regeneration
- › Periodontal tissue in S1 Graft area is esthetic, and it healed stably



Preoperative X-ray



Preoperative Intraoral photo (#12~22 temporary crown)




X-ray after extraction of maxillary left, right lateral incisor




Bone loss area in maxillary right lateral incisor



Molding of S1 bone graft material




X-ray after implant placement and S1 bone graft material use



Bone formation pattern at 4 months after S1 bone graft material use



Final prosthetics



X-ray at 11 months after S1 use



# Open Socket Management in Posterior Area

**Initial assesment** Periodontitis in maxillary right 1st and 2nd premolar and 1st molar areas has progressed for long time with inflammation on the surrounding alveolar bone. There is significant alveolar bone loss on CBCT. Patient wanted implant treatment.


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives


- › Confirmation of alveolar bone loss and remaining bone morphology due to periodontitis in the maxillary right 1st and 2nd premolar and 1st molar area
- › Immediate implant placement after #14, #16 extraction and 3-unit prosthesis production scheduled
- › S1 bone graft material will be used for the bone defect of in the extraction site

## Conclusions

- › Though it was sutured in an open wound state immediately after implant placement with bone graft, secondary healing occurred without significant loss of S1 graft material.
- › In consideration of soft tissue contraction and separation of bone graft material particles when implanted immediately after extraction, transmucosal bone graft that grafts bone graft material up to the top of the soft tissue.
- › Well-formed new bone around the implant was observed on the X-ray at 3 months after surgery.



# Treatment Steps



Inflammation area and alveolar bone shape on preoperative CBCT



Preoperative Intraoral photo



Confirmation of affected bone loss area after extraction of #14, #15, #16



Hydration of S1 bone graft material



Implant placement after S1 graft on broad bone loss area



Implant placement after S1 graft on broad bone loss area



Suture



Intraoral photo at 2 months after surgery



Final prosthetics



X-ray at 3 months after S1 bone graft material use



CBCT at 3 months after S1 bone graft material use

# Narrow Ridge Augmentation in Anterior Area

**Initial assesment** Visited with peri-implantitis detected in another clinic. After removing the existing implant, S1 bone graft material was used for the bone defect.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Implant removal due to peri-implantitis in mandibular anterior (#31)
- › Horizontal bone augmentation surgery performed on the labial side of the residual bone

## Conclusions

- › Final prosthetic placement completed at 3 months after S1 implantation
- › Aesthetic periodontal tissue detected after second implant placement



Preoperative X-ray



Preoperative Intraoral photo



Application of S1 bone graft material for horizontal bone augmentation after implant (1 stage surgery)



Suture



Abutment placement for prosthetics at 3 months after S1 placement



Abutment placement before prosthetics



Gingiva tissue state



Final prosthetics



X-ray after final prosthetics



# Easy Augmentation of Dehiscence Defect

**Initial assesment** 67-year-old woman visited the hospital for reconstruction of a defect in the maxillary left lateral incisor

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (CGF) ☐ No

## Objectives

- › Labial resorption of surrounding alveolar bone after extraction of the maxillary left lateral incisor.
- › Horizontal bone augmentation was planned in consideration of aesthetics, and bone graft material is needed to increase the height of the alveolar bone in the marginal area. With the consent of the patient, blood was collected to process CGF that would be used as membrane.

## Conclusions

- › Despite the lack of bone width in the defect area and in the anatomical structure that made it difficult to maintain the graft material, the desired shape could be easily realized due to the moldable characteristics of S1.
- › The bone graft material position was stably maintained without the use of a membrane, which helps to maintain the shape.



Preoperative Intraoral photo of alveolar bone in the defect area



Incision and flap elevation



Dehiscence defect condition after implant placement



Application of S1 bone graft material



Molding to fit the surrounding alveolar bone shape



Formation of membrane by using CGF



Postoperative X-ray

# Simplified Vertical Augmentation for Advanced Bone Resorption

**Initial assesment** The patient, with extracted 2nd premolar and 1st and 2nd molar at mandibular right posterior region due to severe periodontitis, is planning to receive implant placement along with vertical bone augmentation of the defect area.


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (COLLA) ☐ No

## Objectives

- › The form of the residual bone looks very irregular after removing teeth and the inflammatory granulation tissue due to severe periodontitis, and, accordingly, this situation requires vertical bone augmentation.
- › Induct new bone regeneration by using the S1 bone graft material after implant placement on the #45, #46, #47 regions

## Conclusions


- › Though there is a wide range of bone loss at defect area, it can stably apply S1 bone graft material to the desired region.
- › Without using a non-absorbable shape-fixing membrane, the location of the S1 bone graft material was retained, which helped the regeneration of a new bone.
- › The periodontal tissue of the S1 graft area was placed esthetically and stably.



Preoperative image of the alveolar bone in the defect area



Incision and flap elevation




Alveolar bone defect after implant placement




Application of the S1 bone graft material on the bone defect area



Use of COLLA membrane for vertical bone augmentation



Postoperative X-ray



X-ray at 4 months after bone graft material use



# Dough Staged Sticky Bone Graft

**Initial assesment** A 54-year-old woman visited for the purpose of reconstruction of the mandibular right 2nd premolar and 1st molar defects

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No

## Objectives


- › As the height of the alveolar bone in the mandibular left defect area is not constant, the shape of the marginal area must be reconstructed.
- › Intramarrow penetration is conducted to help the bone graft material create new blood vessels as plenty of time has passed since the extraction, while the remaining bone surface is healed with a dense bone.

## Conclusions

- › On radiographs 3 months after surgery, it was observed that the shape of the marginal bone was harmoniously connected with the adjacent teeth, and the new bone density was also satisfactory.
- › Confirmed that the moldable property of S1 is applied as an advantage in maintaining the 1 wall defect of the distal of the #46 implant placement site.



Postoperative X-ray



Preoperative image of the alveolar bone in the defect area




Incision and flap elevation



Alveolar bone defect after implant placement



Application of S1 on the bone defect area



Formation of membrane by using COLLA



Final coverage using CGF



Suture



X-ray at 3 months after surgery

# Immediate Implant Placement

**Initial assesment** Visited for the reason that the maxillary left central incisor was shaking

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (COLLA) ☐ No

## Objectives

- › Immediate implant placement was performed to preserve the alveolar bone and soft tissue.
- › Labial bone was lost, so after dehiscence defect surgery, the incision was minimized in consideration of edema.


## Conclusions

- › Used S1 and COLLA, the preservation of periodontal tissue was better than before extraction, and esthetically satisfactory results were obtained.
- › The procedure was performed by minimizing the periodontal tissue incision, and the membrane was fixed well, resulting in a good prognosis with solid sutures.

Before



After




Before




After



Before



After






# Treatment Steps



Postoperative X-ray



Preoperative Intraoral photo



Incision and flap elevation



Application of S1 bone graft material after Implant placement



Application of S1 bone graft material



Use of COLLA membrane



Suture



Postoperative X-ray



Intraoral photo at 2 months after surgery



Final prosthetics connection at 3 months after surgery



Intraoral photo after prosthetics



X-ray at 3 months after surgery

# Easy 3D Rigde Augmentation

**Initial assesment** A 75-year-old woman visited for the purpose of restoration of the defect area in the mandibular right 1st premolar and molar


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (CGF) ☐ No

## Objectives

- › A socket-shaped, 3-wall defect was found in the #44 tooth defect after soft tissue elevation during surgery
- › Vertical bone augmentation was performed on #44~46 using S1 bone graft material

## Conclusions

- › S1 bone graft material can be stably applied to the desired area, though the defect area bone width was narrow and the defect area was wide
- › After surgery, it was observed that the shape of the marginal bone was harmoniously connected with the adjacent teeth, and new bone density was also satisfactory.



Alveolar bone shape in the defect area in #44, #45, #46, #47



Implant placement




Application of S1 bone graft material to the bone defect area



Application of S1 bone graft material to the bone defect area



Use of CGF for vertical bone augmentation surgery



Postoperative X-ray



# Labial Fenestration Defect

**Initial assesment** A 57-year-old female patient visited for implant placement of the mandibular right central incisor

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Postoperative X-ray



Preoperative Intraoral photo (#41)



Defect area after implant placement



Application of S1 bone graft material



Application of COLLA membrane



Suture



Postoperative X-ray



Final prosthetics connection at 3 months after surgery



Postoperative Intraoral photo after final prosthetics

# Minimally Invasive Bone Grafting in Anterior Area

**Initial assesment** A 67-year-old female patient who was using partial mandibular anterior dentures visited for the production of a fixed prosthesis. Planned to perform narrowed alveolar bone augmentation by minimizing the incision area during surgery.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Postoperative X-ray



Preoperative Intraoral photo



Application of COLLA membrane after incision (#42)



Application of COLLA membrane after incision (#42)



Application of S1 bone graft material between the membrane and the alveolar bone (#32)



Suture



Intraoral photo after prosthetics



Intraoral photo after prosthetics



X-ray after final prosthetics



# Horizontal Bone Augmentation

**Initial assesment** A 53-year-old female patient visited for the reconstruction of the second premolar and molar area of the mandible. A long period of time has passed after tooth extraction, and residual bone resorption has progressed in the buccal direction. Planned to perform horizontal guided bone regeneration after implant placement

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Postoperative x-ray (#45, #46, #47)



Preoperative Intraoral photo



Implant placement



Application of S1 bone graft material and COLLA membrane



Suture



Intraoral photo after prosthetics



Final prosthetics

# Labial Dehiscence in Upper Anterior Area

**Initial assesment** Visited for restoration of the maxillary left canine loss. Implant placement planned.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



1 Preoperative Intraoral photo



2 Elevation after soft tissue incision



3 Application of 0.15g of S1 bone graft material



4 Molding of S1 bone graft material to fit defect area



5 Application of COLLA membrane



6 Suture



7 Intraoral photo after prosthetics



8 Final prosthetics



# Labial Fenestration in Lower Anterior area

**Initial assesment** Visited for restoration of the missing teeth of the mandibular left central incisor and the right lateral incisor. Planned to perform guided bone regeneration at the same time as implant placement due to the loss of the buccal width of the residual bone in the anterior region.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



1 Preoperative Intraoral photo (#31, #42)



2 Elevation after soft tissue incision



3 Implant placement



4 Application of S1 bone graft material



5 Application of COLLA membrane



6 Suture

# Peri-Implantitis Treatment with laser therapy

**Initial assesment** A 62-year-old female patient, Inflammation next to prosthesis in #37, the implant was placed 10 years ago

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Check the bone defect around the implant in #37 (the left mandibular second molar) in panorama
- › Treatment with laser and filling the bone defect with S1

## Conclusions

- › Removing granuloma on the implant surface and filling in bone defect using S1
- › S1 has been well attached to the implant surface
- › In the future, new bone is expected to be formed well between implant threads.



Preoperative X-ray (#37)



Incision by Laser



Cyst and inflamantion Tissue



Application of S1 bone graft material to the bone defect area



Postoperative x-ray after using S1 bone graft material



# Extraction socket management

**Initial assesment** #27 Patient's back and joint disease due to vertical root fracture (VRF)  
**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

### Objectives

- › After extraction of #27 tooth with vertical root fracture (VRF), immediate implant placement is conducted after removing inflammation, and S1 is applied.

### Conclusions

- › In #27, Alveolar Ridge Preservation (ARP) was performed using S1 without membrane. After the surgery, the soft tissue takes to heal well.
- › After 1 year, it was observed that it was very stable without resorption of marginal bone.



Preoperative (#27)



Socket size measurement after tooth extraction



Immediate implant placement



Hydration of S1



Socket Management by applying S1



Healing abutment placement




#27 Suture



After 1 week



After 1 month & soft tissue healing



After 1 year

# Moldable Augmentaion in Posterior Area

**Initial assesment** Visited for chronic inflammation and tooth mobility of the mandibular left first molar teeth and mandibular right first molar teeth



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extraction plan with inflammation and vertical mobility of both first molar teeth
- › Observation of 3 months after removal of inflammatory tissue after tooth extraction

## Conclusions

- › First molars on both sides of the mandibular were close to the nerve and had a wide range of bone defects, but the stabtivity of new bone regeneration and healing was confirmed
- › Vertical bone augmentation was performed without a membrane, and as a result of CBCT confirmation, stable alveolar bone shape was maintained and new bone regeneration was good.






## Treatment Steps



1  
3 months after extraction of both mandibular first molar teeth



2  
IntraOral photo of left mandibular first molar teeth before surgery



3  
Preoperative CBCT of mandibular right first molar teeth and mandibular left first molar teeth



4  
Mandibular left first molar tooth incision



5  
Mandibular left first molar tooth implantation of implant




6  
Transplantation of S1 bone graft material shaped according to the size of the mandibular left first molar tooth defect



7  
Suture



8  
Postoperative (4 weeks)



9  
Postoperative (4 weeks)



10  
Mandibular left first molar tooth 10 weeks after Healing Abutment tightened and sutured



11  
Mandibular left first molar final prosthesis



12  
IntraOral photo of right mandibular first molar tooth before surgery

## Treatment Steps



Right mandibular first molar gum incision



Right mandibular first molar implant placement



Transplantation of S1 bone graft material shaped according to the size of the defect in the right mandibular first molar




Suture



Post OP (2 months)



2nd surgery



CBCT after bilateral mandibular first molar surgery



Mandibular right first molar prosthesis



Mandibular right first molar prosthesis



# Immediate implant placement

**Initial assesment** Visited for a fractured maxillary central incision that had undergone nerve treatment



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Immediate implant placement after extraction of fractured maxillary central incisor
- › Bone grafting plan with extraction socket and thin labial bone

## Conclusions

- › Bone grafting using Tunneling technique on the thin Labial Bone in the apical region of the maxillary central incisor
- › After 3 months of surgery, even without the membrane, it was observed that the bone was not absorbed and the shape of the bone was maintained well.




# Treatment Steps



Preoperative Panorama



Preoperative intraoral photo



Preoperative CBCT



After the extraction



S1 bone graft



Implant placement in the optimal position



Tunneling technique \_ soft tissue incision



Tunneling technique \_ bone graft material in the apical region




Tunneling technique \_ Suture



Postoperative Panorama



1 month after surgery



CBCT after surgery



## Treatment Steps



Soft tissue healing 2 months after surgery



X-ray picture 2 months after surgery




Soft tissue healing 3 months after surgery



Final Prosthesis Panorama



Maxillary left lateral incisor final prosthesis



Final Prosthesis CBCT

# Immediate implant placement

Dr. Ho Yeul Jang

**Initial assesment** Visited for complaining of discomfort due to inflammation and mobility of the bridge teeth prosthesis of the maxillary central incisor and maxillary lateral incisors



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extraction the maxillary anterior bridge teeth, and remove granulation tissue
- › Planning the implant placement simultaneously with horizontal bone augmentation by observing the absorbed shape on the labial side

## Conclusions

- › Bone grafting and immediate implant placement by removing the granulation tissue of the maxillary central incisors and maxillary lateral incisors
- › After 3 months of surgery, even without the membrane, it was observed that the gum and bone graft shape volume was maintained well






# Treatment Steps



Preoperative Panorama



Preoperative Oral photo



Preoperative CBCT



Extraction



Modeled S1 bone graft material



S1 bone graft



Implant placement in the optimal position




Suture and bone graft with Tunneling technique on the apical part of the upper right lateral incisor



Bone grafting using Tunneling technique on the dentate end of the left maxillary central incisor



Postoperative Panorama



Postoperative CBCT



Immediate loading

## Treatment Steps



Post OP (6 weeks)




Post OP (3 months)



Final Prosthesis



Final Prosthesis Panorama



Final Prosthesis CBCT



# Sinus augmentation

**Initial assesment** Visited for maxillary posterior implant placement

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Maxillary inflammation on the left maxillary first molar and left maxillary second molar requires removal of inflammation
- › Planning for maxillary sinus lift and immediate implant placement by closing the sinus perforation area with PRF

## Conclusions

- › Inflammation of the left maxillary molar is removed with a syringe
- › 3 months after the surgery, the size of inflammation decreased, new bone regeneration and implants were stably maintained.
- › 1 year, 1 year and 6 months after surgery, no inflammatory changes, and observed mature bone tissue.



Preoperative Panorama



Preoperative intraoral photo



Preoperative CBCT



Removal of inflammation using a syringe



Hydration and mixing of S1 powder bone graft material



S1 bone graft material placed on maxillary left 1st molar and maxillary left 2nd molar



Implant placement on maxillary left 1st molar and maxillary left 2nd molar



Healing Abutment Fastening



Postoperative Panorama

## Treatment Steps




Final prosthesis after 3 months




Postoperative CBCT




Final Prosthesis CBCT



Final prosthesis CBCT



Post OP (1 year)



Post OP (1 year and 6 months)



# Socket Management in Posterior Area & Immediate implant placement

**Initial assesment** An 80-year-old man visiting for rapid reconstruction of the food impaction, mobility, smell of inflammation on mandibular right first molar.



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› Immediate implant placement and bone grafting plan after tooth extraction as the the vertical mobility of the mandibular first molar is 3 degrees are exist.

## Conclusions


› After extraction of the mandibular right first molar, inflammation was removed to confirm the thin buccal-lingual cortical bone and large defect.  
› Bone graft was performed after implant placement due to a large defect in extraction site  
› After 3 months of implant placement, the final prosthesis was completed, and shape maintenance, osseointegration, and stable healing of gingival soft tissue were observed.



# Treatment Steps



Preoperative Panorama



Preoperative X-ray



Preoperative intraoral photo



extraction



immediate implant placement



Placement of S1 bone graft material that has moldable characteristics




Insertion of bone graft material



suture




Post OP (4 weeks)



Post OP (4 weeks) X-ray



Post OP (3 months)



Post OP (3 months) X-ray



# Minor bone augmentation (Maxillary sinus lift revision surgery)

**Initial assesment** A 57-year-old woman, needed implant and inflammatory tissue removal because bone regeneration was not successful after the right maxillary sinus lift surgery due to inflammation  
Planning for implant placement after follow-up the result of reoperation for maxillary sinus elevation


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives


- › remove the implant and granulation tissue due to the bone on the right side of the maxilla was not regenerated and alveolar bone destruction because of Maxillary sinus lift surgery inflammation
- › Planning for new bone regeneration and implant placement using sticky bone graft material S1 bone for the large defect in the maxillary sinus.

## Conclusions

- › When removing granulation tissue due to maxillary sinus lift infection, perforation occurred, so closing with a membrane, and S1 bone graft was performed in the large defect.
- › Although it was not an easy case for perforation of the maxillary sinus and large defects, the use of sticky bone made it possible to easily implement the shape of the defect without causing bone graft material particle separation and infection.
- › 6 months after maxillary sinus elevation and implant placement in the large defect, maxillary sinus inflammation was remarkably reduced, and bone level and implant fixation were observed to be stably maintained.



# Treatment Steps



Previous) Maxillary right 1st molar, maxillary right 2nd molar Sinus lift CBCT



implant removal



Confirmation of large defects after maxillary granulation tissue removal



Close perforated site in the maxillary sinus by membrane



S1 Bone preparation




Molding of bone graft materials to fit large defects



Maxillary sinus elevation bone graft revision surgery




Suture



CBCT confirmed that the dispersion of the bone graft material is limited 4 days after surgery



Oral photo showing bone regeneration through gum incision



Post OP (6 months)



implant placement



Check gingival soft tissue healing



Final prosthesis



# Horizontal Alveolar bone Augmentation

**Initial assesment** Visited for a mandibular first molar implantation due to long-term tooth loss

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

- › both sides teeth are tilted due to long-term tooth loss, and horizontal augmentation surgery is required due to absorption of the buccal margin.

## Conclusions

- › S1 bone graft was performed and horizontal alveolar augmentation was performed by membrane, and after 4 months, it was confirmed that the alveolar bone was well formed.
- › Stable restoration of the final prosthesis after implantation due to the sufficient width of the implant through horizontal augmentation



1 Panoramic view at the first visit



2 #36 Incision



3 Implant placement & cortical perforation



4 S1 Bone Graft, membrane



5 S1 Bone Graft, membrane



6 Panoramic view of the first implant surgery



7 Post OP (4 months 2 weeks)



8 Gingival incision to check for bone regeneration



9 Healing abutment after 2nd surgery



10 Suture



11 Mandibular left first molar final prosthesis setting



12 Panoramic view of the final restoration setting

# Treatment Steps



Panoramic view at the first visit



Preoperative Oral photo



Confirmation of bone resorption after incision



Mandibular right second molar implant placement



Mandibular right 1st molar, mandibular right 2nd molar absorbed downward to bone level



Moldable and sticky S1 bone after hydration



Bone Graft to match defect



membrane insertion



Suture



Post OP (1 months)



Post OP (1 months) Panorama



Post OP (5 months)



Signed healing abutment after 2nd OP



Signed healing abutment after 2nd OP



Final prosthesis Panorama



# Horizontal ridge augmentation

**Initial assesment** The right lateral periphery of the maxillary cavity was absorbed so Horizontal bridge augmentation is performed.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › The alveolar bone on the buccal side of the maxillary right first premolar was very thin, so implant placement was performed along with alveolar bone augmentation.

## Conclusions

- › Horizontal ridge augmentation was performed with S1 bone graft material on the absorbed buccal side, and it was confirmed that alveolar bone and volume were created 4 months later.
- › Healing abutment after 2nd surgery



Marginal absorption horizontally on the buccal side of the maxillary right first premolar



Confirmation of marginal areas absorbed by incision and elevation



S1 Bonegraft preparation



Bonegraft on the buccal margin after implant placement on the maxillary right first premolar



Horizontal ridge augmentation



Suture



stich out (+10 day)



intraoral photo (+4 months)



healing cap proceed after 2nd stage surgery



Panorama after 2nd stage surgery

# Horizontal guided bone regeneration

Dr. Goh

Initial assesment

Products                    S1 bone graft material (Powder type, 0.2~1.0mm)   |   Use of membrane ☐ Yes   ☒ No

Objectives

› GBR was performed due to labial upper absorption of the maxillary left central incisor, and implant placement was performed along with surgery.

Conclusions

› Horizontal guided bone regeneration was performed on the absorbed labial side with S1 bone graft material, and it was confirmed that alveolar bone and volume were created 4 months later.



Oral photo of upper left lateral incisor



It was confirmed that labial absorption was confirmed by gingival incision



Confirmation of implant penetration in the labial thin margin



Labial GBR performed after implant placement on the maxillary left lateral incisor



Suture



Suture



Post OP (3 days)



Post OP (3 days)



CBCT (+3 months)



intraoral photo (+4 months)



# Horizontal Alveolar bone Augmentation

**Initial assesment** Visited for implant placement due to loss of maxillary left posterior teeth

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

› Planning the Maxillary sinus elevation and implant placement due to low alveolar bone height in the left posterior region of the maxillary

## Conclusions

› S1 bone graft was performed and horizontal alveolar augmentation was performed by membrane, and after 4 months, it was confirmed that the alveolar bone was well formed.

› Sufficient alveolar bone was formed for implant placement.



1 Gum incision




2 Drilling of upper left posterior teeth



3 Observation of movement of the sinus membrane of the maxillary left second molar during breathing



4 Elevation of the maxillary sinus and insertion of the Boss bone



5 Radiograph after maxillary sinus lift



6 Implant placement



7 Healing Abutment Fastening



8 Suture



9 Post OP

# Narrow Ridge Augmentation in Posterior Area

**Initial assesment** After extracting the existing tooth in narrow buccal bone, using S1 to reinforce the alveolar bone


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

- › After tooth extraction, place S1 in the socket and proceed with open technique using PRF membrane.


## Conclusions

- › Even though proceed with open technique without using a membrane, the soft tissue was well covered, and bone formation was very good after 6 months.





## Treatment Steps



Preoperative x-ray



Confirmation of affected bone loss area after extraction



Application of S1 bone graft material



Bone covering with PRF membrane



Suture



Suture



Post OP




Post OP



Post OP



Post OP



X-ray after S1 bone graft material use

MedPark



NEW BONE  
NEW LIFE



**Tel** +82 2 554 8776 **E-mail** [biz@medpark.net](mailto:biz@medpark.net)  
**[www.medpark.net](http://www.medpark.net)**