Description:

Neomem® FLEXPLUS Resorbable Collagen Membrane is a white, non-friable, resorbable, single layered, conformable collagen membrane matrix manufactured from purified porcine peritoneum. Neomem® FLEXPLUS is sterilized by gamma irradiation and is supplied sterile, non-pyrogenic, and for single use only.

Intended Uses:

Neomem® FLEXPLUS is a bioresorbable, implantable collagen material that is intended for use in oral surgical procedures as a resorbable membrane material for use in:

- simultaneous use of guided bone regeneration (GBR)-membrane and implants;
- augmentation around implants placed in immediate extraction sockets;
- augmentation around implants placed in delayed extraction sockets;
- localized ridge augmentation for later implantation;
- alveolar ridge reconstruction for prosthetic treatment;
- filling of bone defects after root resection, cystectomy or removal of retained teeth;
- guided bone regeneration in dehiscence defects; and
- guided tissue regeneration procedures in periodontal defects.

Contraindications:

Neomem® FLEXPLUS is contraindicated in patients who have:

- acute infections or contaminated wound in the oral cavity;
- a known allergy to collagen; and/or
- clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment, will prevent safe implantation or likely healing.

Administration:

Neomem® FLEXPLUS is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments.

The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect.

Neomem® FLEXPLUS can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to final placement.

Neomem® FLEXPLUS can be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors.

Neomem® FLEXPLUS should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. The mucoperiosteal flap is sutured over the collagen membrane and the wound should be closed completely to avoid accelerated resorption due to membrane exposure. The membrane is expected to be resorbed in approximately 12 to 16 weeks.

Post-Operative Procedures:

Neomem® FLEXPLUS is completely resorbable and should not be removed. Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peride) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated.

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

Probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

Warnings:

Clinicians should use care in screening their patients for any known allergies to collagen. Hypersensitivity reactions have been noted with the use of other products containing collagen; therefore, the possibility exists of developing a local sensitivity response to Neomem® FLEXPLUS.

Precautions:

As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroid therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, indicating a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. Nor has it been evaluated in pregnant women, children and/or in patients with conditions involving extremely severe defects with little periodontium or bone.

Neomem® FLEXPLUS cannot be resterilized or reused. Open, unused Neomem® FLEXPLUS must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

Do not use if the product sterilization barrier or its packaging is compromised.

Adverse Reactions:

Possible complications that can occur with any dental surgery include infection, swelling of the intraoral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing, resorption or ankylosis, with loss of crestal bone height, pain, or complications associated with the use of anesthesia. Minor discomfort may occur for a few days.

Storage:

The product should be stored at room temperature (15°C/59°F to 30°C/86°F). Avoid excessive heat and humidity.

How Supplied:

Neomem® FLEXPLUS is supplied sterile, non-pyrogenic and for single use only.

One (1) membrane per package

Catalogue Number	Size
NF1520	15 mm x 20 mm
NF2030	20 mm x 30 mm
NF3040	30 mm x 40 mm

Caution:

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or dentist.

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