


# Prelude<sup>one</sup><sup>TM</sup>


## ONE COMPONENT UNIVERSAL ADHESIVE

Prelude One is a light-curable single component self-etch adhesive designed to bond composite resins to most dental surfaces. These surfaces include dentin, enamel, zirconia, alumina, metals, porcelain, lithium disilicate, and composite restorations. As a universal adhesive, Prelude One is truly compatible with light cured, dual cured, and auto cured composites (no activator required).


Shear Bond Strength on Tooth Surface (LC Mode)



Shear Bond Strength on Other Surface (LC Mode)



Shear Bond Strength on Tooth Surface (AC Mode)



Available in 5ml bottle and unitdose

### ADVANTAGES

- Convenient one step application
- For use on all surfaces
- No post-op sensitivity
- Extremely high bond strengths to all direct and indirect restorations
- Compatible with dual and self cure composites without the use of an additional activator
- 51/2 micron film thickness

- SBS to Lithium disilicate (E-Max CAD, Ivoclar) with HF treatment (MPa)
- SBS to precious metal (Panceram 62, Pentron) MPa
- SBS to non-precious metal (Cocr) MPa
- SBS to Zirconia (3M Lava) MPa

**ZD ZEST DENTAL SOLUTIONS<sup>TM</sup>**  
 ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS  
 2875 LOKER AVENUE EAST  
 CARLSBAD, CA 92010  
 USA

TEL: 800.262.2310  
 FAX: 760.743.7975  
 EMAIL: SALES@ZESTDENT.COM  
**WWW.ZESTDENT.COM**